

# Technical note

## J.T. Baker® BAKERBOND® PolyPEI™ multimode weak anion exchange chromatography resin

### INTRODUCTION

Utilizing a unique surface chemistry that combines primary, secondary, and tertiary weak anionic exchange sites, BAKERBOND® PolyPEI™ multimode anion exchange chromatography media provides unique selectivity to enable higher purification performance across a wide range of biopharmaceuticals.

Designed and manufactured by Avantor to the high standards established by our J.T. Baker® brand, the BAKERBOND® PolyPEI™ resin provides better selectivity than conventional weak anion exchange media with equivalent capacities. It is capable of separating proteins and peptides having similar isoelectric points (pI). Its fully porous and crosslinked methacrylate spherical beads provide excellent mechanical and chemical stability and are optimized for maximum capacity, high mass transfer, and relatively low backpressure. The resin is provided in a non-hazardous, non-flammable storage solution through Avantor's global supply chain.

Used with Avantor's proven J.T. Baker® family of process chromatography buffers and additives, the BAKERBOND® PolyPEI™ resin can deliver greater efficiencies and higher purity profiles to biopharma chromatography schemes.

### FEATURES


- Unique selectivity delivered by its proprietary mixed mode functionality
- High ionic capacity to improve separation and enable operation at high salt concentration
- Spherical polymethacrylate beads enable uniform packing and provide mechanical stability and chemical resistance for ease-of-use
- Delivery in a non-hazardous, non-flammable buffer solution to eliminate burdensome handling requirements.



## UNIQUE SELECTIVITY

BAKERBOND® PolyPEI™ anion exchange resin is inherently well suited to discriminate between extremely similar molecules. Its unique surface chemistry is able to discern minute differences in primary structure resulting in a unique selectivity, improved separation performance, and higher efficiency downstream operations. As illustrated in Fig.1, BSAB-Lactoglobulin B, and B-Lactoglobulin A having close isoelectric points (pI 4.7, 5.34, and 5.21) are well resolved and can be easily separated.

Fig.1. Separation of Proteins on PolyPEI



**Column:** 0.77 × 10 cm

**Buffer A:** 50 mM TRIS pH 8.0

**Buffer B:** Buffer A plus 1 M NaCl, pH 8.0


**Linear gradient** from 0 to 100% B in 10 CV

**Flow rate:** 1.2 mL/minute

**Sample volume:** 1.5 mL.

Protein	pI	MW	Sample conditions
Lysozyme	10.0	14.7 kD	0.4 mg/mL
Human IgG	7.0	150.0 kD	4.0 mg/mL
BSA	4.7	66.4 kD	2.5 mg/mL
B-Lactoglobulin B	5.3	18.3 kD	1.0 mg/mL
B-Lactoglobulin A	5.2	18.4 kD	1.0 mg/mL


The ability to discern these minute molecular differences is driven by the resin's multimode functionality. The multimode functionality of PolyPEI™ is obtained by covalently bonding Polyethyleneimine (PEI) to the surface of highly cross linked polymethacrylate beads. Polyethyleneimine contains primary, secondary and tertiary nitrogens, which each function as weak anionic exchange sites.



Tertiary, secondary and primary amine anion exchange groups

## HIGH IONIC CAPACITY

Our BAKERBOND® line of polymeric resins were all designed with higher ionic capacities than the industry leading ion exchange and hydrophobic interaction chromatography (HIC) resins. The high ionic capacity increases the salt gradient needed to separate similar molecules, thus improving selectivity. This also enables operations at high salt concentrations. A comparison of ionic capacity between BAKERBOND® PolyPEI™ and several other leading competitive hydrophobic interaction resins can be seen below.




## EASE OF USE

Our BAKERBOND® line of polymeric resins are based on a highly cross linked rigid spherical polymethacrylate media that provides optimal porosity, high mechanical strength, and excellent chemical resistivity delivered with a consistent narrow particle size distribution. These features enable robust resin application, narrow elution bands for concentrated product fractions, and provide consistency and stability to downstream operations.

As illustrated in Fig.2 below, BAKERBOND® PolyPEI™ has an average pore size of 500 Å for an exclusion limit of  $1 \times 10^6$  dalton. This open porosity coupled with the bead's mechanical strength allows for the use of PolyPEI™ at high linear velocities due to good mass transfer.


Fig.2. SEM image of polymethacrylate spherical beads at 2,000 magnification



The average pore size of the beads is 500 Å, providing an exclusion limit of  $1 \times 10^6$  Dalton

A small average particle size of 35 µm combined with a narrow particle size distribution ( $d_{60}/d_{10}$ ) produces narrow elution bands of highly concentrated product for improved efficiency. This increases resolution and decreases pool volumes compared to media based on larger particles. High efficiency in combination with high selectivity enables higher loading levels while maintaining separation.

Fig.3. Particle size distribution of a typical lot of polymethacrylate beads




A median particle size of 35 micron and a particle size distribution of less than 1.4 ( $d_{60}/d_{10}$ ) show a lack of fines and result in high efficiency at low back pressure.

The hydrophilic backbone of BAKERBOND® PolyPEI™ has low non-specific binding and the ability to withstand prolonged contact with commonly used cleaning and sanitizing solutions. PolyPEI™ media can be easily packed to bed heights of up to 40 cm and operated in the linear velocity range of up to 500 cm/hr using conventional columns with pressure ratings of 2 to 10 bar. Our entire line of BAKERBOND® chromatography resins, buffers, and additive are manufactured under the strictest controls and is supported by world class quality systems and application resources.

## STABILITY AND CONSISTENCY

The performance of BAKERBOND® PolyPEI™ remains consistent across a broad range of storage conditions. As illustrated in Fig.4 and Fig.5, little to no change was identified in protein binding capacity BSA after prolonged storage under both acidic and basic conditions. Product performance was evaluated by measuring the dynamic binding capacity of BSA. Similarly, the purification behavior of PolyPEI™ showed little to no change in its ability to separate a mixture of lysozyme, rabbit gamma globulin, and alpha lactalbumin after subjecting the column to a comparable range of storage conditions

Fig.4. Relative Dynamic Binding Capacity (DBC) of BSA under various storage conditions



Column: 7.75 × 100 mm

Binding buffer: 20 mM Sodium Acetate, pH 6.2

Elution buffer: 1 M Sodium Acetate, pH 6.2.


A: Fresh lot in Sterilant buffer at room temperature

B: in 10mM H3PO4 at room temperature for 4.5 years

C: in 0.1N NaOH at room temperature for 4.5 years

D: in Sterilant buffer at low temperature for 4.5 years

Fig.5. Chemical Stability of PolyPEI



The unique selectivity of PolyPEI is not affected by extreme pH conditions over extended periods of time.

**Column:** 10 × 100 mm, sample injected: 5 mL

**Binding buffer:** 20 mM Sodium Acetate, pH 6.2

**Elution buffer:** Binding buffer plus 1 M NaCl, pH 6.2  
**Sample:** lysozyme, rabbit gamma globulin and alpha lactalbumin

**Red:** PolyPEI for 10 months in storage solution at 4 °C

**Blue:** PolyPEI for 10 months in storage solution at room temperature

**Pink:** PolyPEI for 10 months in 10 mM phosphoric acid

**Green:** PolyPEI in 0.1 M NaOH at room temperature for 10 months

**Cyan:** PolyPEI washed dynamically with 0.5 M NaOH for 48 hours at 1.0 mL/min.

## OPERATIONAL FLEXIBILITY

BAKERBOND® PolyPEI™ scales easily to the production environment with reproducible results in columns ranging from 20 cm to 45 cm in internal diameter. The narrow particle size distribution and mechanical stability of the media maintains conventional low pressure column performance at bed heights up to 20 cm. Modern medium pressure columns, rated at 7 bar, can be used to pack beds up to 40 cm in height and enable operation at high linear velocities.

As illustrated in Fig.6, the linear pressure-flow curve for BAKERBOND® PolyPEI™ indicates that there is no resin compression at a bed height of 21.5 cm and a column diameter of 20 cm. In addition, the media does not need to be defined prior to the initial packing or subsequent packings, reducing the consumption of buffers and simplifying operations.

Fig.6. Pressure-Flow relation of polymethacrylate




Fig.7. Efficiency of column (20 id x 21.5 cm high) packed with polymethacrylate based media





Fig.8. Scale – up to a 45 cm ID Column



## CONVENIENT PACKAGING

BAKERBOND® PolyPEI™ resin comes conveniently packaged in a non-hazardous, non-flammable storage buffer, eliminating many of the burdensome shipping, handling, and storage requirements common among chromatography resins.

### Product ordering information

Size	Product Number
1 mL Lab Columns (5-pack)	6059-07
5 mL Lab Columns (5-pack)	6059-25
50 mL	7585-01
100 mL	7585-02
500 mL	7585-03
1 L	7585-04
5 L	7585-05

## PROCESS AND APPLICATION SUPPORT

Avantor has deep expertise in process chromatography optimization and can work with you to help ensure that BAKERBOND® PolyPEI™ delivers the improved performance you need in downstream processing. Technical support from our scientists and application specialists is available from our multiple global research and innovation centers.

### Product Information Summary

#### Characteristics

Functionality	Weak Anion Exchanger
Functional Group	-CH <sub>2</sub> CH <sub>2</sub> NH <sub>2</sub> , (-CH <sub>2</sub> CH <sub>2</sub> ) <sub>2</sub> NH, (CH <sub>2</sub> CH <sub>2</sub> ) <sub>3</sub> N
Ion Exchange Capacity	0.3 – 0.5 Cl-meq/mL
Median Particle Size	35 µm
Particle Size Distribution (d <sub>60</sub> /d <sub>10</sub> )	< 1.4
Median Pore Size	500 Å
Exclusion Limit	1 × 10 <sup>6</sup> Daltons
Operating pH Range	4.5 – 14.0
Cleaning pH Range	1.0 – 14.0
Chemical Stability	All commonly used aqueous buffers, sodium hydroxide, acetic acid, phosphoric acid, guanidine hydrochloride, up to 100% ethanol, methanol, or 2-Propanol.
Shipping Solvent	Media is shipped as 1:1 slurry in acetate buffer containing benzyl alcohol at pH 4.5 and can be stored for up to 5 years at 4–15 °C. The media can also be stored in 20% ethanol or 0.1 M NaOH for up to 5 years.



As our channel brand, VWR, part of Avantor, offers an integrated, seamless purchasing experience that is optimized for the way you do business so you can obtain what you need to move your processes forward. For more information or to buy our products, please visit your local VWR website.

#### AUSTRIA

VWR International GmbH  
Graumanngasse 7  
1150 Vienna  
Tel.: +43 1 97 002 0  
info.at@vwr.com

#### BELGIUM

VWR International bv  
Researchpark Haasrode 2020  
Geldenaaksebaan 464  
3001 Leuven  
Tel.: +32 (0) 16 385 011  
vwr.be@vwr.com

#### CANADA

VWR International  
2360 Argentia Road  
Mississauga, Ontario L5N 5Z7  
Tel.: +1 800 932 5000  
Canada\_Orders@vwr.com

#### CHINA

VWR (Shanghai) Co., Ltd  
Bld.No.1, No.3728 Jinke Rd,  
Pudong New District  
Shanghai, 201203- China  
Tel.: 400 821 8006  
info\_china@vwr.com

#### CZECH REPUBLIC

VWR International s. r. o.  
Veetee Business Park  
Pražská 442  
CZ - 281 67 Stříbrná Skalice  
Tel.: +420 321 570 321  
info.cz@vwr.com

#### DENMARK

VWR International A/S  
Tobaksvejen 21  
2860 Søborg  
Tel.: +45 43 86 87 88  
info.dk@vwr.com

#### FINLAND

VWR International Oy  
Valimotie 9  
00380 Helsinki  
Tel.: +358 (0) 9 80 45 51  
info.fi@vwr.com

#### FRANCE

VWR International S.A.S.  
Le Périgares – Bâtiment B  
201, rue Carnot  
94126 Fontenay-sous-Bois cedex  
Tel.: 0 825 02 30 30\* (national)  
Tel.: +33 (0) 1 45 14 85 00 (international)  
info.fr@vwr.com  
\* 0,18 € TTC/min + prix appel

#### GERMANY

VWR International GmbH  
Hilpertstraße 20a  
D - 64295 Darmstadt  
Tel.: 0800 702 00 07\* (national)  
Tel.: +49 (0) 6151 3972 0 (international)  
info.de@vwr.com  
\*Freecall

#### HUNGARY

VWR International Kft.  
Simon László u. 4.  
4034 Debrecen  
Tel.: +36 52 521130  
info.hu@vwr.com

#### INDIA

VWR Lab Products Private Limited  
No.139, BDA Industrial Suburb,  
6th Main, Tumkur Road, Peenya Post,  
Bangalore, India – 560058  
Tel.: +91 80 28078400  
vwr\_india@vwr.com

#### IRELAND

VWR International Ltd  
Orion Business Campus  
Northwest Business Park  
Ballycoolin  
Dublin 15  
Tel.: +353 (0) 1 88 22 222  
sales.ie@vwr.com

#### ITALY

VWR International S.r.l.  
Via San Giusto 85  
20153 Milano (MI)  
Tel.: +39 02 3320311  
info.it@vwr.com

#### KOREA

VWR International –  
17, Daehak 4-ro, Yeongtong-gu,  
Suwon-si, Gyeonggi-do  
Tel.: +82 31 645 7256  
saleskorea@avantorsciences.com

#### THE NETHERLANDS

VWR International B.V.  
Postbus 8198  
1005 AD Amsterdam  
Tel.: +31 (0) 20 4808 400  
info.nl@vwr.com

#### MEXICO

VWR International, S.de R.L. de C.V.  
Km. 14.5 Carretera  
Tlalnepantla-Cuautitlán  
Col. Lechería  
Tultitlán Edo. de México  
CP 54940  
Tel.: +52 (55) 5005 0100  
vwrmx@vwr.com

#### MIDDLE EAST & AFRICA

VWR International FZ-LLC  
Office 203, DSP Lab Complex,  
Dubai Science Park,  
Dubai, United Arab Emirates  
Tel.: +971 4 5573271  
Info.mea@vwr.com

#### NORWAY

VWR International AS  
Brynsalleen 4,  
0667 Oslo  
Tel.: +47 22 90 00 00  
info.no@vwr.com

#### POLAND

VWR International Sp. z o.o.  
Limbowa 5  
80-175 Gdansk  
Tel.: +48 58 32 38 200  
info.pl@vwr.com

#### PORTUGAL

VWR International -  
Material de Laboratório, Lda  
Centro Empresarial de Alfragide  
Rua da Indústria, nº 6  
2610-088 Amadora  
Tel.: +351 21 3600 770  
info.pt@vwr.com

#### SINGAPORE

VWR Singapore Pte Ltd  
18 Gul Drive  
Singapore 629468  
Tel.: +65 6505 0760  
sales.sg@vwr.com

#### SPAIN

VWR International Eurolab S.L.U.  
C/ Tecnología 5-17  
A-7 Llinars Park  
08450 - Llinars del Vallès  
Barcelona  
Tel.: +34 902 222 897  
info.es@vwr.com

#### SWEDEN

VWR International AB  
Fagerstagatan 18a  
163 94 Stockholm  
Tel.: +46 (0) 8 621 34 00  
kundservice.se@vwr.com

#### SWITZERLAND

VWR International GmbH  
Lerzenstrasse 16/18  
8953 Dietikon  
Tel.: +41 (0) 44 745 13 13  
info.ch@vwr.com

#### UK

VWR International Ltd  
Customer Service Centre  
Hunter Boulevard - Magna Park  
Lutterworth  
Leicestershire  
LE17 4XN  
Tel.: +44 (0) 800 22 33 44  
uksales@vwr.com

#### UNITED STATES

VWR International, LLC  
100 Matsonford Road  
Building One Suite 200  
Radnor, PA 19087  
Tel.: +1 800 932 5000  
VWRCustomerService@vwr.com